

Premi[®]Scan

Instructions to use

With Premi[®]Scan it is possible to objectively measure positive and negative test results of Premi[®]Test by using a desktop scanner.

Requirements to operate:

- Operating System: Windows98 or better.
- Scanner: Recommended HP Scanjet 7400 series (also applicable HP Scanjet 4300 and 6300 series).

Other:

- 20MB of free hard drive space, 32MB of internal memory (RAM).
- A connection to the internet to download a Kodak Q60 Calibration Target File and to enable links to internet sites. Available at: <ftp://ftp.kodak.com/GASTDS/Q60DATA>.
- Black cloth: delivered with the Premi[®]Scan.
- Kodak Professional Q-60 Color Input Target (Cat 190 7914): delivered with the Premi[®]Scan.
- Templates for the ampoules: delivered with the Premi[®]Scan.

Optional:

- Adobe Acrobat Reader to view the manual (downloadable at www.adobe.com)

Please read the scanner operating instructions before performing any measurements.

User License

DSM PremiTest B.V.

Version: 0310
Date of issue: March 17, 2006

DSM Nutritional Products
DSM Premi[®]Test B.V., P.O. Box 1163,
6160 BD Geleen, The Netherlands, Tel. +31 (0)46 4763573
Fax. +31 (0)46 4763737, Internet site: www.premitest.com

Although the facts and suggestions in this publication are based on our own research and are believed reliable, we cannot assume any responsibility for performance or results obtained through the use of our products herein described, nor do we accept any liability for loss or damages directly or indirectly caused by our products. The user is held to check the quality, safety and all other properties of our product prior to use. Nothing herein is to be taken as permission, inducement or recommendation to practise any patented invention without a license.