



Challenges and Experiences in Using New Materials

Bob Ward

President/CEO, DSM PTG

Joseph Chinn

Owner, J Chinn LLC

Jennifer Neff


CEO/CTO, Allvivo Vascular

Janine Robinson

VP Research, Articulinx, Inc.


Ming Wu

VP Eng., Edwards Lifesciences




New Materials in Orthopedic Implants

Dynamic Spinal Fixation



Cervical Spinal Disc



Lumbar Spinal Disc



Hip Joint

Integration of Material and Device Development

Lab and Pilot-Scale Synthesis

Production-Scale Synthesis




Characterization and QC



Component Fab.

Device Assembly



1. Why use a new material?

- **Satisfy bulk and surface property requirements of a device**
- **Obtain a performance advantage (over competitive devices)**
- **Replace an unreliable or unwilling supplier**
 - **Improve lot-to-lot consistency**
 - **Support a chronic implant application**
- **Improve device manufacturing process**
 - **Reduce COGs**
- **Strengthen IP position**
- **Facilitate regulatory approval**
 - **Establish(ed) Master File**
 - **ISO and cGMP Quality Systems**

How does FDA view new materials?

2:40 – 3:20 P.M. Room 210 B

Risk Management and Assessment in Sourcing and Working with New Materials

Jennifer A. Neff, CEO, CTO Allvivo Vascular, Inc.

Evaluation and Testing Within a Risk Management Process: ISO 10993-1:2009 (E)

Table A.1 — Evaluation tests for consideration

Medical device categorization by			Biological effect							
Category	nature of body contact (see 5.2) Contact	contact duration (see 5.3) A – limited (≤ 24 h) B – prolonged (> 24 h to 30 d) C – permanent (> 30 d)	Cytotoxicity	Sensitization	Irritation or intracutaneous reactivity	Systemic toxicity (acute)	Subchronic toxicity (subacute toxicity)	Genotoxicity	Implantation	Haemocompatibility
Surface device	Skin	A	X ^a	X	X					
		B	X	X	X					
		C	X	X	X					
	Mucosal membrane	A	X	X	X					
		B	X	X	X					
		C	X	X	X		X	X		
	Breached or compromised surface	A	X	X	X					
		B	X	X	X					
		C	X	X	X		X	X		
External communicating device	Blood path, indirect	A	X	X	X	X				X
		B	X	X	X	X				X
		C	X	X	X	X	X	X		X
	Tissue/bone/dentin	A	X	X	X					
		B	X	X	X	X	X	X	X	
		C	X	X	X	X	X	X	X	
	Circulating blood	A	X	X	X	X				X
		B	X	X	X	X	X	X	X	X
		C	X	X	X	X	X	X	X	X
Implant device	Tissue/bone	A	X	X	X					
		B	X	X	X	X	X	X	X	
		C	X	X	X	X	X	X	X	
	Blood	A	X	X	X	X	X		X	X
		B	X	X	X	X	X	X	X	X
		C	X	X	X	X	X	X	X	X

^a The crosses indicate data endpoints that can be necessary for a biological safety evaluation, based on a risk analysis. Where existing data are adequate, additional testing is not required.

2. What can the material supplier do to help?

- **Data and Design Support**
 - Liability
- **Samples**
- **Processing Assistance**
- **Support of Regulatory Process**
 - Quality Systems
 - Material Master File

3. How can you help your material supplier?

- **Would you consider supplying your (device) biocompatibility test results to the material supplier for inclusion in *confidential* FDA Material Master Files?**


4. Evolution or Revolution?

- **Is continual improvement and/or modification of established biomaterials preferable to the use of entirely new materials?**
 - **Surface Modification**
 - Improved biostability
 - **Property Enhancement**
 - **Processing Improvements**
 - Gel Reduction
 - Better Injection Molding
 - **Lot-to-lot Consistency**
 - **Expanding Master File**

Improved Biostability: Fluorocarbon End Groups Reduce Environmental Stress Cracking of otherwise identical 80A Polyether-urethanes in Intramuscular Rabbit Implants @ 400 % Strain



Polyether-urethane Control:
3-month explant



**Polyether-urethane with
fluorocarbon end groups:**
6-month explant


Thank You

**DSM Biomedical
Hall D
Booth #2449**

www.dsmpgtg.com


**DSM PTG
2810 7th Street
Berkeley, CA 94710
(510)841-8800**

Evolution vs Revolution



Evolution of Thermoplastic Biomedical Polyurethanes:

Continual Improvement of Existing
Material Platforms



Continuous Improvement of Thermoplastic Polyurethanes for Chronic Implants

Platform Compositions

- **Aromatic (MDI) Polyether-urethane**
 - Strong and hydrolytically stable
- **Aromatic (MDI) Polycarbonate-urethanes**
 - Very strong and oxidatively stable (when ether free)

DSM PTG Enhancements

Via Composition Changes

- **Mixed Soft Segments: Silicone-urethane copolymers with enhanced biostability (ca. 1989)**
- **Use of surface activity and self assembly (of end groups) for surface modification:**
- **Ether-free soft segments**
 - Improved oxidative stability
 - Increased toughness

Continuous Improvement of Thermoplastic Polyurethanes for Chronic Implants

DSM PTG Enhancements

Via Processing Improvements:

- **Replace batch synthesis with continuous synthesis by reactive extrusion**
 - **On-line feedback for MW control**
 - **Reduced particulates**
 - **Rapid development of custom polymers**
- **Gel reduction technologies for higher yields**
 - **During polymer synthesis**
 - **During extrusion of tubing**

Some DSM PTG Thermoplastic Polyurethanes

Bionate[®] polycarbonate-urethane

- **CarboSil[®]** silicone-polycarbonate-urethane
- **CarboSil[®] AL** aliphatic silicone-polycarbonate-urethane
- **Bionate[®] II** polycarbonate-urethane (PCU) with SAME[®] technology


Elasthane[™] polyether-urethane with low MW wax

- **PurSil[®]** silicone-polyether-urethane
- **PurSil[®] AL** aliphatic silicone-polyether-urethane
- **Elasthane[™] II** polyether-urethane with SAME[™] technology*

Note: All polymer families have FDA Master Files

* Developmental material

38 Years of Applications: Polyurethanes in Blood Pumps



AVCO IAB 1971:
*First Clinical
Cardiac Assist
Device*



**Jarvik III
TAH**



Abiocor TAH



Abiomed BVS




**Thoratec PVAD
and IVAD**


TCS HeartMate® LVAD



**Thoratec
HeartMate II**



Jarvik 2000



**Pacemakers
(Various)**

**Sunshine Heart
C-Pulse™ 2009**




1.a. Improved Passive Thrombo-Resistance: *In Vivo* Platelet Thrombus on Catheters is Reduced Three Orders of Magnitude by a Polymeric Surface Modifying Additive (SMA)

Experiment in Progress



Radiolabelled Platelet Uptake vs time




1000-fold reduction in Adherent Platelets on Modified TPU Catheter

1.b. Improved Passive Thrombo-Resistance: Aherent Thrombi on In Vivo Sac-Type LVAD is Reduced with Surface Modifying Additive in Segmented Polyurethane-urea

VAD Blood Sacs: 30-day Calf Explants



SPU with Surface Modification



SPU without Surface Modification

Three-Steps to Antimicrobial Tubing: *Tubing Extruded from Pre-dried 'SAME Polymer' Pellets Made by Continuous Synthesis*

1. Continuous Synthesis of Pellets




2. Dry Polymer Pellets



Tubing Ready for Use




3. Tubing Extrusion and Collection



Self Assembling End Groups:

Self Assembly on Extruded Polyurethane Tubing, e.g. 'Bionate® II'



DSM PTG
Part of DSM Biomedical

SAME™ Technology

Self Assembling Monolayer End Groups

PolymerTechnology_DSM_SAME_NEW

Not to scale

New Device Possibilities with End-Group-Modified Thermoplastics

- Biomaterials must satisfy bulk *and* surface property requirements
- PUs with *Self Assembling Monolayer End Groups (SAME™)* are a 'Biomaterials Toolkit' for independently optimizing bulk *and* surface properties
- SAME™ technology has the potential to provide a simple and economical method for manufacturing tailored (e.g., antimicrobial) polymer surfaces without negative impact on polymer bulk properties.
- Antimicrobial and other SAME polymers can be easily integrated into medical device manufacturing when used as:
 - Neat Structural Materials
 - Surface-modifying additives to other base polymers
 - Coatings on other materials